



## ACETOACETATE - AUT

### DETERMINATION OF ACETOACETATE IN EDTA BLOOD


#### Enzymatic method

Suitable for all analyzers – 150 tests




Product insert with instructions for automated and manual procedures

Stability reagents > 4 years after production

Acetoacetate controls available



#### Settings for automatic analyzers



**Linearity: 1.5 mmol/l**

**Mean CV's: 1.99%**

**Mean recovery: 100.5%**

Product name	Product no.	Quantity
Acetoacetate AUT Reagent Set	3061	30 -150 tests
Acetoacetate Control - AUT, Medium Level	3062	10 x 1 ml
Acetoacetate Control - AUT, High Level	3063	10 x 1 ml
Acetoacetate Control - AUT, Extra High Level	3064	10 x 1 ml





## ACETOACETATE AUTOMATED – ENZYMATIC METHOD

### DETERMINATION OF ACETOACETATE IN EDTA BLOOD

- Enzymatic method
- Suitable for all analyzers – 150 tests
- Product insert with instructions for automated and manual procedures
- Stability reagents > 4 years after production
- Acetoacetate controls available
- Wavelength 340 nm




Products	Product no.	Quantity
<b>Acetoacetate AUT Reagent Set</b>	3061	30 - 150 tests
<b>Acetoacetate Control - AUT, Medium Level</b>	3062	10 x 1 ml
<b>Acetoacetate Control - AUT, High Level</b>	3063	10 x 1 ml
<b>Acetoacetate Control - AUT, Extra High Level</b>	3064	10 x 1 ml

## SUMMARY

### PRINCIPLE

The procedure utilizes the enzyme, D-3-hydroxybutyrate dehydrogenase, which catalyses the following reaction:



In the presence of excess NADH, substantially all acetoacetate is converted to D-3-hydroxybutyrate. The reduction of absorbance at 340 nm due to oxidation of NADH to NAD<sup>+</sup> becomes a measure of the amount of acetoacetate originally present.

### SAMPLE MATERIAL

Deproteinized EDTA blood. Plasma and serum cannot be used.

### LINEARITY

Up to 1.5 mmol/l

### EXPECTED VALUES

*Fasting venous EDTA blood:*  
< 0.1 mmol/l

### QUALITY CONTROL

Products	Product no.	Quantity
<b>Acetoacetate Control - AUT, Medium Level</b>	3062	10 x 1 ml
<b>Acetoacetate Control - AUT, High Level</b>	3063	10 x 1 ml
<b>Acetoacetate Control - AUT, Extra High Level</b>	3064	10 x 1 ml

### QUANTITY OF DETERMINATIONS

Procedure

- Automated : 150 tests
- Manual : 30 tests

### NOTES

1. For in vitro diagnostic use only.
2. For professional use only.
3. Contact INstru**chemie** for the complete validation report and the latest edition product insert.



## CONCENTRATION MEASUREMENTS

The concentrations of a medium, high and extra high sample were measured with an automatic analyzer in order to verify acceptable absorbances.

### Acetoacetate measurements

	Medium	High	Extra High
<b>Δ Absorbance</b>	0.1771	0.3018	0.4998
<b>Concentration (mmol/l)</b>	0.20	0.65	1.29

## TEST CONDITIONS


All tests were performed under the following conditions:

Temperature	: 37 °C
Wavelength	: 340 nm
Light path	: 0.6 cm
Blank	: Distilled or deionized water
Sample	: Deproteinized EDTA blood

## LINEARITY

The Acetoacetate AUT assay is linear up to 1.5 mmol/l.

### Linearity measurements with an automatic analyzer



## SENSITIVITY

The sensitivity (limit of detection) was determined by measuring deproteinized human control material (Acetoacetate concentration = 0 mmol/l) 20 times.

Sensitivity = 3 x standard deviation = 0.06 mmol/l

## PRECISION

The precision is determined by measuring Acetoacetate Control AUT High Level and a deproteinized human blood sample 10 times a day (repeatability) for 5 consecutive days (reproducibility), using an automatic analyzer.

### Repeatability:

	Sample (mmol/l)	Control (mmol/l)
Mean	0.28	0.24
Standard deviation	0.007	0.005
Variation coefficient (%)	2.50	2.08

### Reproducibility:

	Sample (mmol/l)	Control (mmol/l)
Mean	0.28	0.24
Standard deviation	0.007	0.005
Variation coefficient (%)	2.50	2.08

## RECOVERY

The recovery is determined by measuring the Acetoacetate concentration in spiked deproteinized human blood 15 times using an automatic analyzer.

### Recovery:

Added Acetoacetate (mmol/l)	Measured (mmol/l)	Recovery (%)
0.39	0.40	102.6
0.79	0.80	101.3
1.18	1.15	97.5

Printed in the Netherlands, October 2015, Edition 3.0

**FAX – FORM: +31-(0)596 – 634 755****TELEPHONE NO.: +31-(0)596 – 634 831****E-MAIL: info@instruchemie.nl**

<b>Attn</b>	INstruChemie B.V. Sales & Marketing Fons Nieland, Arie Jan de Graaf
<b>Tel no.</b>	+31-(0)596 – 634 831
<b>Fax no.</b>	+31-(0)596 – 634 755
<b>E-mail</b>	info@instruchemie.nl
<b>Date</b>	February 2015
<b>Reference</b>	2151027-1.JB

**I would like:**

<input type="checkbox"/>	<b>Personal contact.</b>	<b>About:</b>
--------------------------	--------------------------	---------------

<input type="checkbox"/>	<b>More information.</b>	<b>About:</b>
--------------------------	--------------------------	---------------

<input type="checkbox"/>	<b>Actual prices.</b>	<b>Of:</b>
--------------------------	-----------------------	------------

<input type="checkbox"/>	<b>Catalogue.</b>	
--------------------------	-------------------	--

**Send to:**

Name:			
Initials		Titel:	M / F
Function:			
Company/Lab:			
Department:			
Address:			
City:			
Telephone:			
Fax:			
Email:			