



## ALDOLASE

### DETERMINATION OF FRUCTOSE-1.6-BISPHOSPHATE ALDOLASE (EC 4.1.2.13) IN SERUM OR PLASMA

#### Enzymatic method

Suitable for all analyzers – 500 tests


Product insert with instructions for automated and manual procedures

Stability > 4 years after production


Aldolase controls available

**Stability working solution:**  
4 weeks (2 – 6°C)  
2 months (-20°C)

#### Settings for automatic analyzers



#### Manual procedure



#### Precision (inter/intra)

	Repeatability (Inter)	Reproducibility (Intra)
Mean	6.55 U/l	6.55 U/l
Standard deviation	0.18 U/l	0.12 U/l
Variation coefficient	2.75 %	1.83 %

**Linearity: 45 U/l**

**Mean CV's: 1.77%**

**Mean recovery: 98.1%**

**Correlation compared to other manufacturers: 0.977**

Product name	Product no.	Quantity
Aldolase Reagent Set	2904	40 -500 tests
Aldolase Calibrator	3008	1 x 1 ml
Aldolase Control Normal Level	2908	10 x 1 ml
Aldolase Control High Level	3007	10 x 1 ml





## ALDOLASE

### DETERMINATION OF FRUCTOSE-1.6-BISPHOSPHATE ALDOLASE (EC 4.1.2.13) IN SERUM OR PLASMA




- Enzymatic Method
- Suitable for all analyzers – 500 tests
- Product insert with instructions for automated and manual procedures
- Stability > 4 years after production
- Aldolase controls available
- Wavelength 340, 334, 365 nm



Products	Product no.	Quantity
<b>Aldolase Reagent Set</b>	2904	40 manual 500 automated
<b>Aldolase Calibrator</b>	3008	1 x 1 ml
<b>Aldolase Control Normal Level</b>	2908	10 x 1 ml
<b>Aldolase Control High Level</b>	3007	10 x 1 ml

## SUMMARY

### ASSAY PRINCIPLE



### SAMPLE MATERIAL

Serum, heparinized plasma or EDTA plasma

### LINEARITY

Up to 45 U/l

### NORMAL VALUE

Up to 7.6 U/l (37°C)

Up to 3.1 U/l (25°C)

### QUALITY CONTROL

Products	Product no.	Quantity
<b>Aldolase Control Normal Level</b>	2908	10 x 1 ml
<b>Aldolase Control High Level</b>	3007	10 x 1 ml

### QUANTITY OF DETERMINATIONS

According to product insert of Aldolase Reagent Set (2904):

Procedure

- Manual : 40 tests
- Automated : 500 tests

### NOTES


1. For in vitro diagnostic use only.
2. For professional use only.
3. Always contact INstruChemie for the complete product insert and latest edition.



## ACTIVITY MEASUREMENTS

The activities of a normal (1), high (2) and extra high (3) human serum were measured with a spectrophotometer.

### Aldolase kinetic measurements



## TEST CONDITIONS

All tests were conducted under the following conditions:

Temperature : 37 °C  
 Wavelength : 340 nm  
 Light path : Analyzer: 0.7 cm / Manual: 1.0 cm  
 Blank : Reagent blank  
 Sample : Serum

## SENSITIVITY


The sensitivity (limit of detection) was determined by measuring human control material (Aldolase activity = 0 U/l) 20 times.

Sensitivity = 3 x standard deviation = 3 x 0.17 = 0.51 U/l

## LINEARITY

The Aldolase assay is linear up to 45 U/l.

### Linearity measurements with an automatic analyzer



## RECOVERY

The recovery is determined by measuring the aldolase activity in spiked human sera 10 times using an automatic analyzer.

### Recovery:

Added aldolase (U/l)	Measured (U/l)	Recovery (%)
4.12	4.05	98.3
10.41	10.20	98.0
25.94	25.39	97.9

## PRECISION

The precision is determined by measuring a human serum and Aldolase Calibrator 10 times a day (repeatability) for 5 consecutive days (reproducibility), using an automatic analyzer.

### Repeatability:

	Sample (U/l)	Calibrator (U/l)
Mean	6.55	12.47
Standard deviation	0.18	0.17
Variation coefficient	2.75	1.36


### Reproducibility:

	Sample (U/l)	Calibrator (U/l)
Mean	6.55	12.48
Standard deviation	0.12	0.14
Variation coefficient	1.83	1.12

## CORRELATION

Pearsons' correlation is determined by measuring the Aldolase activity in multiple human sera with reagent of INstruChemie and reagent from another manufacturer.

### Correlation measured with an automatic analyzer (U/l)



**FAX – FORM: +31-(0)596 – 634 755****TELEPHONE NO.: +31-(0)596 – 634 831****EMAIL: info@instruchemie.nl**

<b>Attn</b>	INstruchemie B.V. Sales & Marketing Fons Nieland, Arie-Jan de Graaf
<b>Tel no.</b>	+31-(0)596 – 634 831
<b>Fax no.</b>	+31-(0)596 – 634 755
<b>E-mail</b>	info@instruchemie.nl
<b>Date</b>	February 2015
<b>Reference</b>	2150205-2.FEN

**I would like:**

<input type="checkbox"/>	<b>Personal contact.</b>	<b>About:</b>
--------------------------	--------------------------	---------------

<input type="checkbox"/>	<b>More information.</b>	<b>About:</b>
--------------------------	--------------------------	---------------

<input type="checkbox"/>	<b>Actual prices.</b>	<b>Of:</b>
--------------------------	-----------------------	------------

<input type="checkbox"/>	<b>Catalogue.</b>	
--------------------------	-------------------	--

**Send to:**

Name:			
Initials		Title:	M / F
Function:			
Company/Lab:			
Departement:			
Address:			
City:			
Telephone:			
Fax:			
Email:			