



β-HBA - AUTOMATIC

DETERMINATION OF D-3-HYDROXYBUTYRATE (β-HBA) IN EDTA BLOOD

Enzymatic method

Suitable for all analyzers – 150 tests




Product insert with instructions for automated and manual procedures

Stability reagents > 5 years after production

β-HBA/Lactate/Pyruvate controls available

Acetoacetate, Lactate,
Pyruvate and β-HBA
from 1 sample

Settings for automatic analyzers



Linearity: 18 mmol/l

Mean CV's: 1.74 %

Mean recovery: 101.6 %

Product name	Product no.	Quantity
β-HBA AUT Reagent Set	3076	30 -150 tests
β-HBA/Lactate/Pyruvate Control - AUT, Low-Normal Level	3112	10 x 1 ml
β-HBA/Lactate/Pyruvate Control - AUT, High Level	3113	10 x 1 ml
β-HBA/Lactate/Pyruvate Control - AUT, Extra High Level	3114	10 x 1 ml





β-HBA AUTOMATED – ENZYMATIC METHOD

DETERMINATION OF D-3-HYDROXYBUTYRATE (β-HBA) IN EDTA BLOOD

- Enzymatic method
- Suitable for all analyzers – 150 tests
- Product insert with instructions for automated and manual procedures
- Stability reagents > 5 years after production
- β-HBA/Lactate/Pyruvate controls available
- Wavelength 340 nm




Products	Product no.	Quantity
β-HBA - AUT Reagent Set	3076	30 -150 tests
β-HBA/Lactate/Pyruvate Control - AUT, Low-Normal Level	3112	10 x 1 ml
β-HBA/Lactate/Pyruvate Control - AUT, High Level	3113	10 x 1 ml
β-HBA/Lactate/Pyruvate Control - AUT, Extra High Level	3114	10 x 1 ml

SUMMARY

PRINCIPLE

The procedure utilizes the enzyme, D-3-hydroxybutyrate dehydrogenase, which catalyses the following reaction:



In the presence of excess NAD^+ , substantially all β-hydroxybutyrate (β-HBA) is converted to acetoacetate. The increase in absorbance at 340 nm, due to the reduction of NAD^+ to NADH, is directly proportional to the β-HBA concentration in the sample.

SAMPLE MATERIAL

Deproteinized EDTA blood. Plasma and serum cannot be used.

LINEARITY

Up to 18 mmol/l

EXPECTED VALUES

Fasting venous EDTA blood:

Adults: 0 – 0.43 mmol/l

Children: 0.02 – 0.30 mmol/l

QUALITY CONTROL

Products	Product no.	Quantity
β-HBA/Lactate/Pyruvate Control - AUT, Low-Normal Level	3112	10 x 1 ml
β-HBA/Lactate/Pyruvate Control - AUT, High Level	3113	10 x 1 ml
β-HBA/Lactate/Pyruvate Control - AUT, Extra High Level	3114	10 x 1 ml

QUANTITY OF DETERMINATIONS

Procedure

- Automated : 150 tests
- Manual : 30 tests

NOTES

1. For in vitro diagnostic use only.
2. For professional use only.
3. Contact INstru**chemie** for the complete validation report and the latest edition product insert.



CONCENTRATION MEASUREMENTS

The concentrations of a normal, high and extra high sample were measured with an automatic analyzer in order to verify acceptable absorbances.


β-HBA measurements

	Normal	High	Extra High
Δ Absorbance	0.0585	0.2456	0.4386
Concentration (mmol/l)	0.50	2.10	3.75

LINEARITY

The β-HBA AUT assay is linear up to 18 mmol/l.

Linearity measurements with an automatic analyzer



PRECISION

The precision is determined by measuring deproteinized human blood sample and β-HBA Control AUT High Level 10 times a day (repeatability) for 5 consecutive days (reproducibility), using an automatic analyzer.

Repeatability:

	Sample (mmol/l)	Control (mmol/l)
Mean	0.31	0.47
Standard deviation	0.005	0.006
Variation coefficient (%)	1.61	1.28

Reproducibility:

	Sample (mmol/l)	Control (mmol/l)
Mean	0.31	0.47
Standard deviation	0.006	0.010
Variation coefficient (%)	1.93	2.13

TEST CONDITIONS

All tests were performed under the following conditions:

Temperature	: 37 °C
Wavelength	: 340 nm
Light path	: 0.7 cm
Blank	: Distilled or deionized water
Sample	: Deproteinized EDTA blood

SENSITIVITY

The sensitivity (limit of detection) was determined by measuring deproteinized human control material (β-HBA concentration = 0 mmol/l) 20 times.

Sensitivity = 3 x standard deviation = 0.01 mmol/l

RECOVERY

The recovery is determined by measuring the β-HBA concentration in spiked deproteinized Bovine Serum Albumin 10 times using an automatic analyzer.


Recovery:

Added β-HBA (mmol/l)	Measured (mmol/l)	Recovery (%)
0.57	0.57	100.0
3.40	3.59	105.6
11.33	11.24	99.2

CORRELATION

Pearsons' correlation is determined by measuring the β-HBA concentration in multiple deproteinized human blood samples with INstruChemie β-HBA AUT (3076) and reagent from another manufacturer

Correlation measured with an automatic analyzer (mmol/l)





FAX – FORM: +31-(0)596 – 634 755

TELEPHONE NO.: +31-(0)596 – 634 831

EMAIL: info@instruchemie**.nl**

Attn	INstru chemie B.V. Sales & Marketing Fons Nieland, Arie Jan de Graaf
Tel no.	+31-(0)596 – 634 831
Fax no.	+31-(0)596 – 634 755
E-mail	info@instru chemie .nl
Date	February 2015
Reference	2150217-2.FEN

I would like:

<input type="checkbox"/>	Personal contact.	About:
--------------------------	--------------------------	---------------

<input type="checkbox"/>	More information.	About:
--------------------------	--------------------------	---------------

<input type="checkbox"/>	Actual prices.	Of:
--------------------------	-----------------------	------------

<input type="checkbox"/>	Catalogue.	
--------------------------	-------------------	--

Send to:

Name:			
Initials		Titel:	M / F
Function:			
Company/Lab:			
Department:			
Address:			
City:			
Telephone:			
Fax:			
Email:			