



FREE FATTY ACIDS - NEFA

DETERMINATION OF FREE FATTY ACIDS (NON ESTERIFIED FATTY ACIDS) IN SERUM OR PLASMA

Enzymatic method

Suitable for all analyzers – 100 tests


Product insert with instructions for automated and manual procedures

Stability > 3 years after production


FFA controls available

Small volumes per vial
=
Minimized excess reagents

Settings for automatic analyzers



Manual procedure



Added FFA (mmol/l)	Measured (mmol/l)	Recovery (%)
1.11	1.11	100.0
2.22	2.23	100.5
3.33	3.34	100.3

Linearity: 4.00 mmol/l

Mean CV's: 2.04%

Mean recovery: 100.3%

Correlation compared to another manufacturer: 0.997

Product name	Product no.	Quantity
Free Fatty Acids Reagent Set	3055	20 -100 tests
Free Fatty Acids Control High Level	3052	10 x 1 ml
Free Fatty Acids Control Normal Level	3053	10 x 1 ml
Free Fatty Acids Control Extra High Level	3054	10 x 1 ml



FREE FATTY ACIDS – NEFA
DETERMINATION OF FREE FATTY ACIDS (NON ESTERIFIED FATTY ACIDS) IN SERUM OR PLASMA


- Enzymatic method
- Suitable for all analyzers – 100 tests
- Product insert with instructions for automated and manual procedures
- Stability > 3 years after production
- FFA controls available
- Wavelength 546 nm



Products	Product no.	Quantity
Free Fatty Acids Reagent Set	3055	20 - 100 tests
Free Fatty Acids Control , High Level	3052	10 x 1 ml
Free Fatty Acids Control , Normal Level	3053	10 x 1 ml
Free Fatty Acids Control , Extra High Level	3054	10 x 1 ml

SUMMARY
PRINCIPLE

Free Fatty Acids react with Coenzyme A in the presence of Acyl CoA Synthetase to Acyl CoA. Acyl CoA is oxidized to 2,3-trans-Enoyl-CoA and H₂O₂. The intensity of the purple pigment is proportional to the concentration of Free Fatty Acids.


SAMPLE MATERIAL

Serum or EDTA Plasma

Stability of the samples: 2 days at 2-6 °C
Store samples at -20 °C when immediate measurement is not possible

LINEARITY

Up to 4.00 mmol/l

EXPECTED VALUES

Serum or EDTA Plasma:

Men : 0.1 – 0.60 mmol/l

Women : 0.1 – 0.45 mmol/l

QUALITY CONTROL

Products	Product no.	Quantity
Free Fatty Acids Control , High Level	3052	10 x 1 ml
Free Fatty Acids Control , Normal Level	3053	10 x 1 ml
Free Fatty Acids Control , Extra High Level	3054	10 x 1 ml

QUANTITY OF DETERMINATIONS

Procedure

- Automated : 100 tests
- Manual : 20 tests

NOTES

1. For in vitro diagnostic use only.
2. For professional use only.
3. Contact INstruChemie for the complete design output report and the latest edition product insert.

CONCENTRATION MEASUREMENTS

The concentrations of a normal high and extra high sample were measured with an automatic analyzer in order to verify acceptable absorbances.


FFA measurements

	Low	Normal	High
Absorbance	0.0258	0.109	0.216
Concentration (mmol/l)	0.19	0.80	1.57

LINEARITY

The FFA assay is at least linear up to 4.00 mmol/l.

FFA linearity measurements with an automatic analyzer



PRECISION

The precision is determined by measuring a serum sample and Control Extra High Level 10 times a day (repeatability) for 5 consecutive days (reproducibility), using an automatic analyzer.

Repeatability:

	Sample (mmol/l)	Control (mmol/l)
Mean	0.24	1.51
Standard deviation	0.006	0.027
Variation coefficient (%)	2.50	1.79

Reproducibility:

	Sample (mmol/l)	Control (mmol/l)
Mean	0.24	1.50
Standard deviation	0.005	0.027
Variation coefficient (%)	2.08	1.80

TEST CONDITIONS

All tests were performed under the following conditions:

Temperature	: 37 °C
Wavelength	: 546 nm
Light path	: 0.7 cm
Blank	: Distilled or deionized water
Sample	: Serum / EDTA plasma

SENSITIVITY

The sensitivity (limit of detection) was determined by measuring basic calibrator material (FFA concentration = 0 mmol/l) 20 times.

$$\text{Sensitivity} = 3 \times \text{standard deviation} = 3 \times 0.003 = 0.009 \text{ mmol/l}$$

RECOVERY

The recovery is determined by measuring the FFA concentration in a spiked human serum sample 10 times using an automatic analyzer.


Recovery:

Added FFA (mmol/l)	Measured (mmol/l)	Recovery (%)
1.11	1.11	100.0
2.22	2.23	100.5
3.33	3.34	100.3

CORRELATION

Pearsons' correlation is determined by measuring the FFA concentration in multiple human sera with INstruChemie reagent and reagent of another manufacturer.

FFA correlation measured with an automatic analyzer (mmol/l)





FAX – FORM: +31-(0)596 – 634 755

TELEPHONE NO.: +31-(0)596 – 634 831

E-MAIL: info@instruchemie**.nl**

Attn	INstru chemie B.V. Sales & Marketing
Tel no.	+31-(0)596 – 634 831
Fax no.	+31-(0)596 – 634 755
E-mail	info@instru chemie .nl
Date	August 2018
Reference	2180816-1.FEN

I would like:

<input type="checkbox"/>	Personal contact.	About:
--------------------------	--------------------------	---------------

<input type="checkbox"/>	More information.	About:
--------------------------	--------------------------	---------------

<input type="checkbox"/>	Actual prices.	Of:
--------------------------	-----------------------	------------

<input type="checkbox"/>	Catalogue.	
--------------------------	-------------------	--

Send to:

Name:			
Initials		Titel:	M / F
Function:			
Company/Lab:			
Department:			
Address:			
City:			
Telephone:			
Fax:			
E-mail:			