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
Products	Product no.	Quantity
Homocysteine Reagent Set	2487	50 – 70 Determinations
Homocysteine Calibration Standard	2478	1 x 2 ml
Homocystein Internal Standard	2492	1 x 2 ml
Homocysteine Control Plasma Normal Level	2485	1 x 2 ml
Homocysteine Control Plasma High Level	2486	1 x 2 ml

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Visiting address  
 Zwet 26  
 9932 AB Delfzijl  
 The Netherlands

Postal address  
 P.O. Box 120  
 9930 AC Delfzijl  
 The Netherlands

**DETERMINATION OF HOMOCYSTEINE BY HPLC-FD IN SERUM OR PLASMA****CLINICAL BACKGROUND AND ASSAY PRINCIPLE**

Homocysteine is a risk factor for premature cardiovascular disease. Determination of homocysteine in plasma and serum is useful both in the diagnosis and follow-up of folate and cobalamin (Vitamin B12) deficiencies and rare inborn errors causing hyper homocysteinemia. The methionine loading test can help diagnosis (2). (Oral 100 mg/kg body weight). Homocysteine is bound to plasma proteins by disulfide bonds. Therefore it is necessary for the determination of total Homocysteine to reduce these disulfide bonds. After deproteinization and bringing the solution to the right pH, derivatization take place.

A good separation of Cysteine, Cysteinylglycine and Homocysteine will be obtained with an Inertsil Analytical column Polaris C18-A. Collected blood samples should be centrifuged immediately after collection and the serum or plasma should be stored at -20 °C when not directly analyzed. Before analysis frozen samples must be thawed at room temperature.

**SAMPLING**

Serum or Plasma

**ANALYTICAL CONDITIONS:**

Column : Polaris C18A, 5 µm packing  
Flow Rate : 0.8 ml/min  
Temperature : 25 °C  
Detection : Fluorescence (excitation 390 nm, emission 500 nm)  
Injection volume : 20 µl  
Runtime : Samples 12 min.

**INTERPRETATION**

Peak No.	Component	Polaris C18-A, 5 µm Approx. retention time
1.	Internal Standard	2 min.
2.	Cysteine	4 min.
3.	Cysteinylglycine	5 min.
4.	Homocysteine	7 min.

**LINEARITY**

Up to 150 µmol/l

**EXPECTED VALUES**

	Male	Female	After methionine loading test
Homocysteine, plasma	6.18 – 13.37 µmol/l	3.27 – 12.43 µmol/l	
Homocysteine, serum	8.56 – 17.12 µmol/l	3.92 – 16.84 µmol/l	< 60 µmol/l

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