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
Products	Product no.	Quantity
Homocysteine Reagent Set	2487	50 – 70 Determinations
Homocysteine Calibration Standard	2478	1 x 2 ml
Homocystein Internal Standard	2492	1 x 2 ml
Homocysteine Control Plasma Normal Level	2485	1 x 2 ml
Homocysteine Control Plasma High Level	2486	1 x 2 ml

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Visiting address
 Zwet 26
 9932 AB Delfzijl
 The Netherlands

Postal address
 P.O. Box 120
 9930 AC Delfzijl
 The Netherlands

DETERMINATION OF HOMOCYSTEINE BY HPLC-FD IN SERUM OR PLASMA**CLINICAL BACKGROUND AND ASSAY PRINCIPLE**

Homocysteine is a risk factor for premature cardiovascular disease. Determination of homocysteine in plasma and serum is useful both in the diagnosis and follow-up of folate and cobalamin (Vitamin B12) deficiencies and rare inborn errors causing hyper homocysteinemia. The methionine loading test can help diagnosis (2). (Oral 100 mg/kg body weight). Homocysteine is bound to plasma proteins by disulfide bonds. Therefore it is necessary for the determination of total Homocysteine to reduce these disulfide bonds. After deproteinization and bringing the solution to the right pH, derivatization take place.

A good separation of Cysteine, Cysteinylglycine and Homocysteine will be obtained with an Inertsil Analytical column Polaris C18-A. Collected blood samples should be centrifuged immediately after collection and the serum or plasma should be stored at -20 °C when not directly analyzed. Before analysis frozen samples must be thawed at room temperature.

SAMPLING

Serum or Plasma

ANALYTICAL CONDITIONS:

Column : Polaris C18A, 5 µm packing
Flow Rate : 0.8 ml/min
Temperature : 25 °C
Detection : Fluorescence (excitation 390 nm, emission 500 nm)
Injection volume : 20 µl
Runtime : Samples 12 min.

INTERPRETATION

Peak No.	Component	Polaris C18-A, 5 µm Approx. retention time
1.	Internal Standard	2 min.
2.	Cysteine	4 min.
3.	Cysteinylglycine	5 min.
4.	Homocysteine	7 min.

LINEARITY

Up to 150 µmol/l

EXPECTED VALUES

	Male	Female	After methionine loading test
Homocysteine, plasma	6.18 – 13.37 µmol/l	3.27 – 12.43 µmol/l	
Homocysteine, serum	8.56 – 17.12 µmol/l	3.92 – 16.84 µmol/l	< 60 µmol/l

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