



ZINC

Colorimetric method (5-Br-PAPS)

Zinc standard included


60 - 150 determinations per set

Serum, plasma, cerebrospinal fluid or urine

Serum and urine controls available

With and without deproteinization


Settings for automatic analyzers



Precision (inter/intra)

	Repeatability (Inter)	Reproducibility (Intra)
Mean	18.3 µmol/l	18.3 µmol/l
Standard deviation	0.21 µmol/l	0.19 µmol/l
Variation coefficient	1.15 %	1.04 %

Manual procedure



Linearity: 155 µmol/l

Mean CV's: 1.63%

Mean recovery: 97.8%

Correlation compared to other manufacturers: 0.981

Product name	Product no.	Quantity
Zinc Reagent Set	2299	60 -150 tests
Zinc Serum Control Low Level	3031	1 x 3 ml
Zinc Serum Control High Level	3032	1 x 3 ml
Zinc Urine Control Low Level	3033	1 x 3 ml
Zinc Urine Control High Level	3034	1 x 3 ml





ZINC

DETERMINATION OF ZINC IN SERUM, PLASMA OR CEREBROSPINAL FLUID

- Colorimetric method (5-Br-PAPS)
- For Manual and/or Automated Procedures
- Instrument Application Sheets Available
- Use Serum, Plasma, Cerebrospinal fluid or Urine
- Incl. Zinc Standard
- Wavelength 560 nm



Products	Product no.	Quantity
Zinc Reagent Set	2299	60-150 tests
Zinc Serum Control, Low Level	3031	1 x 3 ml
Zinc Serum Control, High Level	3032	1 x 3 ml
Zinc Urine Control, Low Level	3033	1 x 3 ml
Zinc Urine Control, High Level	3034	1 x 3 ml

SUMMARY

PRINCIPLE

Zinc, in a pH 8.60 buffer system, forms with specific complexant 5-Br-PAPS a stable colored complex. The color intensity of which is proportional to the amount of zinc in the sample. The interferences, due to oligoelements present in the sample, are eliminated using particular reaction condition and specific masking agents.

SAMPLE MATERIAL

Serum, plasma, cerebrospinal fluid and urine. EDTA-plasma and hemolyzed samples cannot be used.

LINEARITY

Up to 155 µmol/l

EXPECTED VALUES

Serum, plasma:

- Male: 10.7 - 17.5 µmol/l
- Female: 11.1 - 19.5 µmol/l

Urine: 52.2 - 122.3 µmol/24 hours

QUALITY CONTROL

Products	Product no.	Quantity
Zinc Serum Control, Low Level	3031	1 x 3 ml
Zinc Serum Control, High Level	3032	1 x 3 ml
Zinc Urine Control, Low Level	3033	1 x 3 ml
Zinc Urine Control, High Level	3034	1 x 3 ml

QUANTITY OF DETERMINATIONS

Procedure

- Manual : 60 tests
- Automated : 150 tests

NOTES

1. For in vitro diagnostic use only.
2. For professional use only.
3. Always contact INstruChemie for the complete product insert and latest edition.



CONCENTRATION MEASUREMENT

The concentrations of a low, normal and high sample were measured with a spectrophotometer in order to verify acceptable absorbance.


Zinc measurements

	Low	Normal	High
Extinction	0.049	0.086	0.255
Concentration (µmol/l)	9.1	15.8	24.4

LINEARITY

The Zinc assay is linear up to 155 µmol/l.

Linearity measurements with an automatic analyzer



PRECISION

The precision is determined by measuring Zinc Serum Control High Level and a human serum 10 times a day (repeatability) for 5 consecutive days (reproducibility), using an automatic analyzer.

Repeatability:

	Sample (µmol/l)	Control (µmol/l)
Mean	18.3	24.8
Standard deviation	0.21	0.53
Variation coefficient (%)	1.15	2.14

Reproducibility:

	Sample (µmol/l)	Control (µmol/l)
Mean	18.3	24.8
Standard deviation	0.19	0.54
Variation coefficient (%)	1.04	2.17

TEST CONDITIONS

All tests were conducted under the following conditions:

Temperature	: 37 °C
Wavelength	: Analyzer: 545 nm / Manual: 560 nm
Light path	: Analyzer: 0.7 cm / Manual: 1.0 cm
Blank	: Reagent blank
Sample	: Serum

SENSITIVITY

The sensitivity (limit of detection) was determined by measuring human control material (Zinc Concentration = 0 µmol/l) 20 times.

$$\text{Sensitivity} = 3 \times \text{standard deviation} = 3 \times 0.2 = 0.6 \mu\text{mol/l}$$

RECOVERY

The recovery is determined by measuring the Zinc concentration in spiked human sera 10 times using an automatic analyzer.


Recovery:

Added Zinc (µmol/l)	Measured (µmol/l)	Recovery (%)
6.8	6.7	98.5
17.0	16.9	99.4
101.7	97.2	95.6

CORRELATION

Pearsons' correlation is determined by measuring the Zinc concentration in multiple human sera with reagent of INstruChemie and reagent from another manufacturer.

Correlation measured with an automatic analyzer (µmol/l)





instru**chemie**[®]

Science without f(r)iction

Zwet 26 9932AB Delfzijl
The Netherlands

Tel. +31 (0)596 - 634 831
Fax +31 (0)596 - 634 755



FAX – FORM: +31-(0)596 – 634 755

TELEPHONE NO.: +31-(0)596 – 634 831

EMAIL: info@instruchemie**.nl**

Attn	INstru chemie B.V. Sales & Marketing Fons Nieland, Arie Jan de Graaf
Tel no.	+31-(0)596 – 634 831
Fax no.	+31-(0)596 – 634 755
E-mail	info@instru chemie .nl
Date	June 2013
Reference	2130612-1.FEN

I would like:

<input type="checkbox"/>	Personal contact.	About:
--------------------------	--------------------------	---------------

<input type="checkbox"/>	More information.	About:
--------------------------	--------------------------	---------------

<input type="checkbox"/>	Actual prices.	Of:
--------------------------	-----------------------	------------

<input type="checkbox"/>	Catalogue.	
--------------------------	-------------------	--

Send to:

Name:			
Initials		Titel:	M / F
Function:			
Company/Lab:			
Departement:			
Address:			
City:			
Telephone:			
Fax:			
Email:			