

**25-OH VITAMIN D3 - CONTROLS**

AS CONTROL FOR ACCURACY AND PRECISION FOR THE DETERMINATION OF 25-OH VITAMIN D3.

Independent Controls
For Spectrophotometric, HPLC, LCMS
Expiration date > 6 years
Lyophilised
3x levels

Value concentration in:
µg/L & nmol/L

SUMMARY**COMPOSITION**

Lyophilised preparations with analytical results for 25-OH Vitamin D3 by HPLC.

APPLICATION

The 25-OH Vitamin D3 Controls are used as an assayed quality control plasma to monitor the accuracy and precision for the determination of 25-OH Vitamin D3 by all methods.

ANALYTICAL RESULTS

Component	Method	Mean Value	Range
25-OH Vitamin D3 Control, Low Level	HPLC	9 µg/l	7 – 11
		24 nmol/l	19 – 29
25-OH Vitamin D3 Control, Normal Level	HPLC	21 µg/l	17 – 25
		53 nmol/l	43 – 63
25-OH Vitamin D3 Control, High Level	HPLC	75 µg/l	65 – 85
		188 nmol/l	163 – 213

PRODUCTS

Products	Product no.	Quantity
25-OH Vitamin D3 Control Low Level	2869	10 x 3 ml
25-OH Vitamin D3 Control Normal Level	2870	10 x 3 ml
25-OH Vitamin D3 Control High Level	2871	10 x 3 ml

**NOTES**

1. For in vitro diagnostic use only.
2. For professional use only.
3. Always contact INstruChemie for the complete product insert and latest edition.

Printed in the Netherlands, August 2018, Edition 4.0