

LIPIDS CALIBRATOR

LYOPHILISED HUMAN SERUM PREPARATION WITH ANALYTICAL RESULTS FOR SEVERAL LIPIDS

Independent Calibrator
 Cholesterol, HDL Cholesterol,
 Triglycerides
 Human Serum Preparation
 Expiration date > 5 years
 Lyophilised



SUMMARY

COMPOSITION

Lyophilised human serum with certified analytical results for Total Cholesterol, HDL-Cholesterol and Triglycerides. By a special preservation procedure the lipoproteins remain native as appeared from density gradient ultracentrifugation and agarose gel electrophoresis.

APPLICATION

As calibrator for Total Cholesterol, HDL-cholesterol and Triglycerides determinations.

ANALYTICAL RESULTS

Product name:

Lipids Calibrator, Medium – High Level – Cholcal

Component:

Total Cholesterol
 HDL-Cholesterol
 Triglycerides

Concentration:

6.6 mmol/l
 1.72 mmol/l
 1.99 mmol/l

PRODUCTS

Products	Product no.	Quantity
Lipids Calibrator, Medium – High Level – Cholcal	2298	10 x 2 ml



NOTES

1. For in vitro diagnostic use only.
2. For professional use only.
3. Always contact INstru chemie for the complete product insert and latest edition.
4. Printed in the Netherlands, July 2015 – version 3.0