



WATER HPLC GRADE

WATER HPLC GRADE

Ready for use
 Expiration date > 2 years
 Store at room temperature




SUMMARY

WATER HPLC GRADE

REF 2195 1 liter 15 - 25 °C
 LOT 50010915 2017-09-30 IVD

non-volatile matter: 0.0001 % non-volatile matter on ignition: 0.0003 %
 conductivity: $0,5 \times 10^{-4} \text{ Ohm}^{-1} \text{ M}^{-1}$ absorbance of largest eluted peak: 0.002 AU

40 ml of water through a C18 reverse phase column and eluting with a linear water / acetonitrile gradient at 5 % per minute at a flow rate of 2 ml per minute



instruchemie****[®]
 Zwet 26 9932 AB Delfzijl The Netherlands

PRODUCTS

Products	Product no.	Quantity
WATER HPLC GRADE	2195	1 x 1 liter



NOTES

1. For in vitro diagnostic use only.
2. For professional use only.
3. Contact INstru**chemie** for the latest edition product insert.

Printed in the Netherlands, August 2018, version 5.0