



## WATER HPLC GRADE

### WATER HPLC GRADE

Ready for use  
 Expiration date > 2 years  
 Store at room temperature




### SUMMARY

**WATER HPLC GRADE**

REF 2195      1 liter      15 - 25 °C  
 LOT 50010915      2017-09-30      IVD

non-volatile matter: 0.0001 %    non-volatile matter on ignition: 0.0003 %  
 conductivity:  $0,5 \times 10^{-4} \text{ Ohm}^{-1} \text{ M}^{-1}$     absorbance of largest eluted peak: 0.002 AU

40 ml of water through a C18 reverse phase column and eluting with a linear water / acetonitrile gradient at 5 % per minute at a flow rate of 2 ml per minute



**instru**chemie****<sup>®</sup>  
 Zwet 26 9932 AB Delfzijl The Netherlands

### PRODUCTS

Products	Product no.	Quantity
WATER HPLC GRADE	2195	1 x 1 liter



### NOTES

1. For in vitro diagnostic use only.
2. For professional use only.
3. Contact INstru**chemie** for the latest edition product insert.
4. Printed in the Netherlands, Oktober 2015, version 4.0